

SURGICAL LENSES

VITRECTOMY | SINGLE USE | ACCESSORIES

VITRECTOMY FOR A RETINAL DETACHMENT WITH MULTIPLE TEARS

A 25-gauge Vitrectomy was performed using the **Miniquad XL SSV**. After clearing the core and peripheral vitreous, PFCL heavy liquid was injected to flatten and stabilize the posterior pole. Vitreous attached to the margins of breaks was also removed. Air fluid exchange was done and endodrainage was carried out through an existing break. Once the retina was flat, 360° endolaser was carried out around the breaks.

Manish Nagpal, MD FRCS FASRS
Director of Retina Foundation,
Ahmedabad, India

INDIRECT VITRECTOMY LENSES


Volk offers a suite of vitrectomy lenses over a range of optical profiles, designed to cater to the full spectrum of vitrectomy procedures with the highest quality Volk optics for the best surgical visualization.

| LENS | FIELD OF VIEW | IMAGE MAG | CONTACT DIAMETER | PRIMARY APPLICATION |
|-----------------|---------------|-----------|------------------------|--|
| HRX Vit Lens | 130° / 150° | 0.43x | 11.35 mm / SSV 16.0 mm | Far-Peripheral Indirect Vitreoretinal Procedures |
| Mini Quad® XL | 112° / 134° | 0.39x | 11.35 mm / SSV 16.0 mm | Indirect Viewing and Treatment of Peripheral Retinal Disorders |
| Mini Quad® | 106° / 127° | 0.39x | 11.35 mm / SSV 16.0 mm | Indirect Viewing and Treatment of Peripheral Retinal Disorders |
| DynaView | 95° / 127° | 0.39x | 8.08 mm | Treatment of Retinopathy of Prematurity |
| Central Retinal | 73° / 88° | 0.71x | 11.35 mm / SSV 16.0 mm | High Magnification Indirect Viewing and Treatment of the Central Retina |
| Super Macula® | 64° / 77° | 1.03x | 11.35 mm | Highest Magnification Indirect Viewing and Treatment of the Central Retina |

HRX Vit Lens




- PRIMARY APPLICATION**
Far-Peripheral Indirect Vitreoretinal Procedures
- + High index glass delivers widest field, distortion-free retinal views of any surgical lens
 - + Small profile ring facilitates instrument manipulation and surgical procedures
 - + Available in standard and patented self-stabilizing contact (SSV) options for best ergonomics
 - + Ideal for retinal detachments, PVR, giant retinal tears and works seamlessly in fluid and air filled eyes
 - + Available in autoclave sterilizable design (see page 52)



Mini Quad® XL



- PRIMARY APPLICATION**
Indirect Viewing and Treatment of Peripheral Retinal Disorders
- + Wide field of view of the entire retina including the ora serrata
 - + Ideal for retinal detachments, giant retinal tears, PDR, including diabetic cases requiring endolaser to the periphery
 - + Available in standard and self-stabilizing contact (SSV) options for best ergonomics



“CRYSTAL CLEAR VISIBILITY & STABILITY

The Volk HRX and MiniQuad XL are my absolute go-to lenses for all my vitrectomy procedures. The wide-field view offered by these lenses allows for crystal clear visibility through all mediums such as fluid, air, PFCL, or silicon oil. Vitrectomy is all about *The View* and these contact lenses provide the best possible view to operate and to get optimum, distortion-free video footage for teaching and academics. Complex cases such as Retinal Detachments with PVR, Giant Retinal Tears, and Diabetic Tractional Detachments have become easier to manage as the Mini Quad XL and HRX lenses provide a seamless view of the extreme periphery to do a thorough clean-up and flatten the retina effectively. The self-stabilizing (SSV) component adds superb stability to this


lens and I don't need any ring or assistant to support it for me. The only time I shift to another lens is when I want to do fine work on the macula like epiretinal membrane peeling or ILM peeling, which is when I move to the Volk Flat SSV lens for that part of the procedure to get the best magnified stereoscopic view of the macula.”

- Manish Nagpal, MD FRCS FASRS
Director of Retina Foundation, Ahmedabad, India

Mini Quad®




- PRIMARY APPLICATION**
Indirect Viewing and Treatment of Peripheral Retinal Disorders
- + Wide field of view of the entire retina including the ora serrata
 - + Smaller ring facilitates manipulation within the orbit
 - + Ideal for retinal detachments, PDR and giant retinal tears
 - + Available in standard and self-stabilizing contact (SSV) options
 - + Available in autoclave sterilizable design (see page 52)



DynaView




- PRIMARY APPLICATION**
Treatment of Retinopathy of Prematurity
- + Enhanced design provides wide field imaging out to the ora serrata
 - + Minified housing facilitates extension of instruments
 - + Reduced contact size ideal for pediatric examination and treatment such as bilateral retinal detachment, vitreous hemorrhage, ROP



Central Retinal




- PRIMARY APPLICATION**
High Magnification Indirect Viewing and Treatment of the Central Retina
- + High resolution, high magnification imaging to the equator
 - + Ideal for epiretinal membranes, diabetic membranes, vitreo macular traction, macular holes, submacular surgeries, and other small detail procedures of the central retina
 - + Available in standard and self-stabilizing contact (SSV) options
 - + Available in autoclave sterilizable design (see page 52)



Super Macula®



- PRIMARY APPLICATION**
Highest Magnification Indirect Viewing and Treatment of the Central Retina
- + High resolution, highest magnification imaging of the central retina
 - + Provides excellent magnification for fine peeling of epiretinal membrane as well as ILM. Ideal for macular holes, vitreo macular traction, and submacula surgeries
 - + 2x field of view compared to plano/concave direct image lenses



AUTOCLAVABLE INDIRECT VITRECTOMY LENSES

| LENS | FIELD OF VIEW | IMAGE MAG | CONTACT DIAMETER | PRIMARY APPLICATION |
|----------------------|---------------|-----------|------------------------|--|
| HRX ACS® | 130° / 150° | 0.43x | 11.38 mm / SSV 16.0 mm | Widest Field Views for Vitreoretinal Procedures |
| Mini Quad® ACS® | 106° / 127° | 0.48x | 11.38 mm / SSV 16.0 mm | Peripheral Indirect Vitreoretinal Procedures |
| Central Retinal ACS® | 73° / 88° | 0.71x | 11.38 mm / SSV 16.0 mm | High Magnification Indirect Vitreoretinal Procedures |


HRX ACS®



PRIMARY APPLICATION

Widest Field Views for Vitreoretinal Procedures

- + Superior high-index glass design ensures widest field views of any vitrectomy lens
- + Advanced aspheric design provides unmatched high resolution imaging
- + Ideal for retinal detachments, PDR and giant retinal tears
- + Steam sterilizable for reduced processing time




Mini Quad® ACS®



PRIMARY APPLICATION

Peripheral Indirect Vitreoretinal Procedures

- + Steam sterilizable for reduced processing time
- + Smaller ring facilitates manipulation within the orbit
- + Ideal for retinal detachments, PDR and giant retinal tears




Central Retinal ACS®



PRIMARY APPLICATION

High Magnification Indirect Vitreoretinal Procedures

- + High resolution, high magnification imaging to the equator
- + Steam sterilizable for reduced processing time
- + Ideal for epiretinal membranes, diabetic membranes, vitreo macular traction, macular holes, submacular surgeries, and other small detail procedures of the central retina



AUTOCLAVABLE SURGICAL BIO LENSES

Combine the optical excellence of Volk lenses with the comfort of reduced processing time in a surgical environment with the autoclavable lens line.

| LENS | FIELD OF VIEW | IMAGE MAG | LASER SPOT MAG | WORKING DISTANCE | RING DIAMETER | PRIMARY APPLICATION |
|----------|---------------|-----------|----------------|------------------|---------------|---|
| 20D ACS® | 46° / 60° | 3.13x | 0.32x | 50 mm | 55.4 mm | Industry Standard Diagnostic Lens in an Autoclavable Format |
| 28D ACS® | 53° / 69° | 2.27x | 0.44x | 33 mm | 45.9 mm | Fundus Scanning Lens in an Autoclavable Format |


20D ACS®



PRIMARY APPLICATION

Industry Standard Diagnostic Lens in an Autoclavable Format

- + Steam sterilizable for use in a surgical environment
- + High quality Permaview™ glass withstands the rigors of repeated sterilization
- + Perfectly corrected for field curvature, astigmatism, and aberrations




28D ACS®



PRIMARY APPLICATION

Fundus Scanning Lens in an Autoclavable Format

- + Steam sterilizable for use in a surgical environment
- + High quality Permaview™ glass withstands the rigors of repeated sterilization
- + Excellent for small pupil diagnosis and treatment including LIO (Laser Indirect Ophthalmoscope)




HIGH RESOLUTION (HR) DIRECT VITRECTOMY LENSES

Volk's High Resolution Direct Image lenses utilize a high-index glass to deliver superior image quality. This robust glass type is highly resistant to the rigors of continued steam sterilization and will not deteriorate or discolor.

Volk's No Stabilizing Ring (NSR) range of lenses allow suitable stability without the need for suturing or stabilizing rings. Two of the lenses in the group are also available in a no suture ring design. The profiles of these two lenses allows them to stabilize suitably without the need for an additional stabilizing ring.

| LENS | FIELD OF VIEW | IMAGE MAG | CONTACT DIAMETER | PRIMARY APPLICATION |
|----------------------|--|---|------------------|--|
| HR Direct 1x | 30° | 1.0x | 11.2 mm | Direct Image Vitreoretinal Surgery of the Central Retina |
| HR Direct Bi-Concave | 45° (Mid Field, Fluid) 30° (AFX, Air) | 0.49x (Mid Field, Fluid) 1.0x (AFX, Air) | 11.2 mm | Wide Field and AFX (Air Fluid Exchange) Direct Image Vitreoretinal Surgery |
| HR Direct High Mag | 20° | 1.35x | 11.2 mm | High Magnification Direct Image Vitreoretinal Surgery of the Central Retina |
| HR Direct 20° Prism | 40° (Offset 20°) | 0.53x | 11.2 mm | Off Axis Wide Field Direct Image Vitreoretinal Surgery |



A case of sub ILM blood collection in which the ILM was peeled to expose the blood, followed by aspiration. The blood is partly whitish in color due to de-hemoglobinization which occurs over time. A Flat SSV Lens was used for this procedure.

- Image courtesy of Dr. Manish Nagpal, Ahmedabad, India


HR Direct 1x



Stabilizing Ring:
VHRD1XACS
No Stabilizing Ring:
VHRD1XNSRACS

PRIMARY APPLICATION Direct Image Vitreoretinal Surgery of the Central Retina

- + High-index glass delivers highest resolution direct image of the central retina
- + Highly suited for repeated steam sterilization with no material degradation
- + Standard design fits all major suture rings
- + Unique optional no stabilizing ring (NSR) design available
- + Ideal for visualizing the posterior pole in ILM peeling





HR Direct Bi-Concave



VHRDBCACS

PRIMARY APPLICATION Wide Field and AFX (Air Fluid Exchange) Direct Image Vitreoretinal Surgery

- + High-index glass in a bi-concave design delivers highest resolution imaging for wide field and AFX procedures
- + Ideal for visualizing fundus through an air filled cavity
- + Highly suited for repeated steam sterilization with no material degradation
- + Standard design fits all major suture rings




HR Direct High Mag



Stabilizing Ring:
VHRDHMACS
No Stabilizing Ring:
VHRDHMNSRACS

PRIMARY APPLICATION High Magnification Direct Image Vitreoretinal Surgery of the Central Retina

- + High-index glass delivers highest resolution, high magnification of the central retina
- + Best suited for detailed work of the macula
- + Highly suited for repeated steam sterilization with no material degradation
- + Standard design fits all major suture rings
- + Unique optional no stabilizing ring (NSR) design available




HR Direct 20° Prism



VHRD20PACS

PRIMARY APPLICATION Off Axis Wide Field Direct Image Vitreoretinal Surgery

- + High-index glass delivers highest resolution off axis (20°) direct image retinal views
- + Improved design delivers wider field (40°) off axis views
- + Highly suited for repeated steam sterilization with no material degradation
- + Ideal for visualizing the posterior peripheral fundus through direct imaging



DIRECT VITRECTOMY LENSES

SELF STABILIZING (SSV)

Volk's Surgical Vitrectomy lenses were developed in collaboration with Dr. K.V Chalam and are available in 7 designs to meet all the visualization needs of a retina surgeon. The SSV® (self-stabilizing) contact element eliminates the need for sutures or rings and provides excellent stability. The compact lens design provides greater spatial access without interfering with instruments.

| LENS | FIELD OF VIEW | IMAGE MAG | CONTACT DIAMETER | PRIMARY APPLICATION |
|---|------------------|-----------|------------------|---|
| Direct Image Flat SSV® (ACS®) | 30° | 0.92x | 11.9 mm | Routine Direct Image Vitreoretinal Surgery of the Central Retina |
| Direct Image High Mag SSV® (ACS®) | 28° | 1.50x | 11.9 mm | High Magnification Direct Image Vitreoretinal Surgery of the Central Retina |
| Direct Image Mid Field SSV® (ACS®) | 40° | 0.50x | 8.0 mm | Wide field of view for pan retinal examination and laser treatments |
| Direct Image 15° Prism SSV® (ACS®) | 30° (15° Offset) | 0.90x | 11.9 mm | Off Axis Direct Image Vitreoretinal Surgery |
| Direct Image 30° Prism SSV® (ACS®) | 30° (30° Offset) | 0.90x | 10.0 mm | Off Axis Direct Image Vitreoretinal Surgery |
| Direct Image 45° Prism SSV® (ACS®) | 30° (45° Offset) | 0.90x | 10.0 mm | Off Axis Direct Image Vitreoretinal Surgery |
| Direct Image AFX SSV® (ACS®) (Air Fluid Exchange - Air Filled Eye) | 30° | 0.82x | 11.9 mm | Direct Image Vitreoretinal Surgery During Air Fluid Exchange Procedures |


Direct Image Flat SSV® ACS®



VFLATSSVACS

PRIMARY APPLICATION
Routine Direct Image Vitreoretinal Surgery of the Central Retina

- + Delivers high resolution direct image of the central retina
- + Steam sterilizable for reduced processing time
- + Most popular lens for high resolution macula work such as epiretinal membrane peeling and ILM peeling




Direct Image High Mag SSV® ACS®



VFHMSSVACS

PRIMARY APPLICATION
High Magnification Direct Image Vitreoretinal Surgery of the Central Retina

- + Delivers high resolution, high magnification direct image of the central retina
- + Steam sterilizable for reduced processing time
- + Ideal for detailed work of the macula with high magnification like macular holes, membrane peeling, tractional retinal detachments




Direct Image Mid Field SSV® ACS®




VMFSSVACS

PRIMARY APPLICATION
Wide Field Direct Image Vitreoretinal Surgery

- + Bi-concave design provides widest field available in a direct image lens
- + Can be used for air/gas exchange procedures
- + Steam sterilizable for reduced processing time




Direct Image 15° Prism SSV® ACS®



VPRISMSSVACS

PRIMARY APPLICATION
Off Axis Direct Image Vitreoretinal Surgery

- + Design delivers 15° off axis retinal views
- + Steam sterilizable for reduced processing time
- + Ideal for direct visualization of the mid-peripheral fundus




Direct Image 30° Prism SSV® ACS®



V30PRISMSSVACS

PRIMARY APPLICATION
Off Axis Direct Image Vitreoretinal Surgery

- + Design delivers 30° off axis retinal views
- + Steam sterilizable for reduced processing time
- + Ideal for direct visualization of the posterior peripheral fundus




Direct Image 45° Prism SSV® ACS®



V45PRISMSSVACS

PRIMARY APPLICATION
Off Axis Direct Image Vitreoretinal Surgery

- + Design delivers 45° off axis retinal views
- + Steam sterilizable for reduced processing time
- + Ideal for direct visualization of the posterior peripheral fundus




Direct Image AFX SSV® ACS®



VAFXSSVACS

PRIMARY APPLICATION
Direct Image Vitreoretinal Surgery During Air Fluid Exchange Procedures

- + Delivers high resolution central retinal imaging
- + Steam sterilizable for reduced processing time
- + Ideal for Air Fluid exchange procedures



SINGLE-USE
SURGICAL BIO
LENSES



Volk®1 Single-Use Surgical BIO lenses combine high-quality optics that Volk is known for and the convenience of pre-sterilization into a ready-to-use design. Volk’s single-use surgical BIO lenses enable convenient pre- and post-operative examination and laser treatment.

These single-use lenses minimize the risk of infection and cross-contamination and reduce the cost and time associate with reprocessing.

Single-use lenses are pre-sterilized and individually-packaged in a Tyvek® pouch. Single-use lenses are sold in boxes of 10.


Volk®1
Single-Use 20D



V20LCD10

PRIMARY APPLICATION
Industry Standard Diagnostic Lens
in a Single-Use Format

- + Perfectly balanced magnification and field of view make this lens ideal for general diagnostic examination
- + Provides excellent views of the optic disc and macula
- + Anti-reflective coating greatly reduces distracting glare




Volk®1
Single-Use 28D



V28LCD10

PRIMARY APPLICATION
Fundus Scanning Lens in a
Single-Use Format

- + Excellent for wide field examination and treatment through a small pupil
- + Compatible with LIO (Laser Indirect Ophthalmoscope)
- + Excellent lens for ROP rounds to reduce infection risk in high-risk babies



“SAFE & EFFICIENT

Since the reports of using reusable lenses during Retinopathy of Prematurity (ROP) screening rounds were linked to infection transmission and serious adverse outcomes in the NICU, I have explored different options to maintain sterile equipment for use during my ROP screening rounds. I feel that the quality and field-of-view of the Volk Single-Use 28D lens is equivalent to the standard and I currently use a separate Volk Single-Use 28D lens for each infant during ROP screening rounds to reduce the risk of infection transmission between infants being examined. I have found that using Volk Single-Use 28D lenses for ROP screening rounds is more efficient than following a protocol to disinfect and reuse standard lenses between infants being screened.”

- S. Grace Prakalapakorn, MD, MPH
Pediatric Ophthalmologist, Durham, NC, USA

SINGLE-USE
DIRECT VITRECTOMY
LENSES



Available in six popular designs, these lenses deliver high resolution direct-image retinal views for all vitrectomy procedures. Designed in collaboration with Dr. K.V. Chalam, the SSV® (self-stabilizing) contact design element eliminates the need for sutures or rings. They are packaged individually in easy-to-open single-use Tyvek® pouches and are boxed in quantities of 10 lenses. These single-use lenses minimize the risk of infection and cross-contamination and reduce the cost and time associate with reprocessing.

| LENS | FIELD OF VIEW | IMAGE MAG | CONTACT DIAMETER | PRIMARY APPLICATION |
|---------------------------------|------------------|-----------|------------------|---|
| Volk®1 Single-Use Flat Standard | 36° | 1.0x | 14.8 mm | Routine Direct Image Vitreoretinal Surgery of the Central Retina |
| Volk®1 Single-Use Flat SSV® | 30° | 0.92x | 16.0 mm | Routine Direct Image Vitreoretinal Surgery of the Central Retina |
| Volk®1 Single-Use Magnifying | 30° | 1.50x | 14.8 mm | High Magnification Direct Image Vitreoretinal Surgery of the Central Retina |
| Volk®1 Single-Use Wide Field | 48° | 0.50x | 14.8 mm | Wide Field Direct Image Vitreoretinal Surgery |
| Volk®1 Single-Use Bi-Concave | 25° | 0.80x | 14.8 mm | Direct Image Vitreoretinal Surgery During Air Fluid Exchange |
| Volk®1 Single-Use 30° Prism | 33° (Offset 30°) | 1.0x | 14.8 mm | Off Axis Direct Image Vitreoretinal Surgery |


Volk®1 Single-Use
Flat Standard



VFD10

PRIMARY APPLICATION
Routine Direct Image Vitreoretinal Surgery
of the Central Retina

- + Ideal for visualizing the posterior pole in ILM peeling
- + Silicone ring base




Volk®1 Single-Use
Flat Self Stabilizing
SSV®



VFLATSSVD10

PRIMARY APPLICATION
Routine Direct Image Vitreoretinal Surgery
of the Central Retina

- + Ideal for visualizing the posterior pole in ILM peeling
- + Patented SSV (self-stabilizing) feet for maximum stability and greater access for instrumentation when working closer to the center axis.




Volk®1 Single-Use Magnifying



VMD10

PRIMARY APPLICATION
High Magnification Direct Image Vitreoretinal Surgery of the Central Retina

- + Ideal for detailed macular work due to high 1.50x magnification
- + Silicone ring base




Volk® Single-Use Wide Field



VWFD10

PRIMARY APPLICATION
Wide Field Direct Image Vitreoretinal Surgery

- + Ideal for wide field imaging of the posterior pole
- + Silicone ring base




Volk®1 Single-Use Bi-Concave



VBCD10

PRIMARY APPLICATION
Direct Image Vitreoretinal Surgery During Air Fluid Exchange

- + Ideal for air-fluid exchange procedures
- + Silicone ring base




Volk®1 Single-Use 30° Prism



V30PD10

PRIMARY APPLICATION
Off Axis Direct Image Vitreoretinal Surgery

- + Ideal for direct visualization of the posterior peripheral fundus
- + Silicone ring base



SURGICAL ACCESSORIES

Suture Ring



VSR52

PRIMARY APPLICATION
Provides a Stable Lens Platform During Vitreoretinal Surgery

- + Premium surgical implant grade titanium for optimal durability and ease of sterilization
- + Larger radius provides enhanced functionality and safety during use
- + Compatible with all Volk direct and indirect contact vitrectomy lenses (except SSV® styles)

Infusion Handle



VINFHAN

PRIMARY APPLICATION
Infusion of Saline Solution Beneath the Lens During Vitreoretinal Surgery

- + Flushes blood and debris providing a clear view during surgery
- + Autoclave sterilizable for reduced processing time
- + Ideal for diabetic surgery

VitreoLens Handle



DynaView Vit, Mini Quad Vit: VVITHAN-LG

Central Retinal Vit, HRX Vit, Super Macula Vit, Mini Quad XL Vit, Central Retinal ACS®, HRX ACS®, Mini Quad® ACS®: VVITHAN-MQXL

PRIMARY APPLICATION
Holding and Stabilization of Lenses During Vitreoretinal Surgery

- + Holds vitrectomy lenses stably to assist during vitreoretinal surgery
- + Malleability allows user to bend the handle to suit their preference
- + Autoclave sterilizable for reduced processing time

Sterilization Tray



Large Tray: VSCB

Small Tray: VS CA

PRIMARY APPLICATION
Sterilization of Ophthalmic Lenses

- + Autoclave safe and approved for use with ETO
- + Small tray (2.7" x 1.5" x 1.25") houses Volk surgical and smaller indirect and slit lamp lenses
- + Large tray (6" x 2.5" x 1.25") houses the largest Volk lenses and accessories including vitrectomy handles

SURGICAL GONIO LENSES

Volk's Surgical Gonioprism lenses leverage the same proprietary optical design and manufacturing principles as Volk's diagnostic lenses. Each surgical gonio lens is designed and tested in partnership with numerous surgeons resulting in the best optical clarity, maximum visualization, surgeon & microscope friendly ergonomics, and optimized for patient comfort.

| LENS | IMAGE MAG | CONTACT DIAMETER | HANDLE LENGTH | PRIMARY APPLICATION |
|---------------------|-----------|------------------|---------------|---|
| VVG Lens | 1.20x | 10.2 mm | 84 mm | Direct Views for Micro-Invasive Glaucoma Surgery (MIGS) and all Intraoperative Gonio Procedures |
| Surgical Gonio Lens | 1.20x | 10.3 mm | 75 mm | Direct Views for Intraoperative Gonio Procedures |

Volk Vold Gonio (VVG) Lens




VTSVVG

PRIMARY APPLICATION
Direct Views for Micro-Invasive Glaucoma Surgery (MIGS) and all Intraoperative Gonio Procedures

- + Thornton-style stabilization ring provides maximum control of the globe
- + Floating ring design minimizes corneal pressure to prevent anterior chamber distortion
- + Visualizes angle in primary phaco position with minimal microscope and head adjustments
- + Designed in collaboration with Dr. Steven Vold and refined with doctors across the world to ensure maximum usability
- + Sterilizable by either steam autoclave or ethylene oxide (ETO)

1.20x
IMAGE MAG



Surgical Gonio Lens




VSGACS

PRIMARY APPLICATION
Direct Views for Intraoperative Gonio Procedures

- + Lightweight titanium handle and chip resistant lens design with adjustable lens orientation
- + Enables clear visualization of the angle for surgery
- + Lens design enables comfortable positioning against the cornea
- + Lens position can be adjusted relative to the handle: for left hand and right hand or center position allowing freedom of movement
- + Applicable for MIGS procedures
- + Sterilizable by either steam autoclave or ethylene oxide (ETO)

1.20x
IMAGE MAG



“SUPERB VISUALIZATION

The Volk surgical gonioprism allows superb visualization of the angle and conforms well to the cornea with minimal coupling agent. The handle is well sized to fit under the increasing size of microscope stacks and the ability to rotate the lens allows additional surgical freedoms while maintaining positional comfort.”

- J. Morgan Micheletti, MD *Cataract, Refractive, & Anterior Segment Surgeon
Berkeley Eye Center, Houston, Texas, USA*

VOLK VOLD GONIO LENS

MICRO-INVASIVE GLAUCOMA SURGERY

A Revolution in MIGS

For maximum control, clearer angle image, and minimal corneal pressure, choose the Volk VVG Lens for Micro-Invasive Glaucoma Surgery (MIGS) and other intraoperative surgical gonio procedures.



Floating Lens

Multiple Degrees of Freedom

Stabilizing Ring

- Stabilize and Control the Globe**
with Thornton-style fixation ring
- Eliminate Anterior Chamber Distortion**
floating lens minimizes pressure on the cornea
- Minimal Microscope & Head Adjustments**
visualize angle in primary phaco position
- Withstands Repeat Sterilization**
compatible with both steam and gas sterilization

SPECIFICATIONS

| IMAGE MAG | CONTACT DIAMETER | RING DIAMETER | HANDLE LENGTH |
|-----------|------------------|---------------|---------------|
| 1.20x | 10.2 mm | 15.2 mm | 84 mm |



“STABILITY FOR MIGS

The floating lens and stabilizing Thornton Ring assist you with rotating the eye so you can easily visualize the trabecular meshwork... and stabilize for perfect visualization.”

- Michael S. Berlin, MD
Director of Glaucoma Institute of Beverly Hills, West Hollywood, CA, USA