




CODICE:	VTSVVG
CND:	Z12120199
RDM:	1886209

DIRECT VIEWS FOR MICRO-INVASIVE **GLAUCOMA SURGERY (MIGS) AND** ALL INTRAOPERATIVE GONIO **PROCEDURES** 

+ THORNTON-STYLE STABILIZATION RING PROVIDES MAXIMUM CONTROL OF THE GLOBE + FLOATING RING DESIGN MINIMIZES CORNEAL PRESSURE TO PREVENT ANTERIOR CHAMBER DISTORTION + VISUALIZES ANGLE IN PRIMARY PHACO POSITION WITH MINIMAL MICROSCOPE AND HEAD ADJUSTMENTS + DESIGNED IN COLLABORATION WITH DR. STEVEN VOLD AND REFINED WITH DOCTORS ACROSS THE WORLD TO ENSURE MAXIMUM USABILITY + STERILIZABLE BY EITHER STEAM AUTOCLAVE OR ETHYLENE OXIDE (ETO)

1.20x IMAGE



LENS	IMAGE MAG	CONTACT DIAMETER	HANDLE LENGTH	PRIMARY APPLICATION
VVG Lens	1.20x	10.2 mm	84 mm	Direct Views for Micro-Invasive Glaucoma Surgery (MIGS) and all Intraoperative Gonio Procedures
Surgical Gonio Lens	1.20x	10.3 mm	75 mm	Direct Views for Intraoperative Gonio Procedures



WWW.SIROFTALMICA.COM